

```
// Generated: Wednesday, June 04, 2014 10:00:30 AM
```

```
profile MSK
```

```
{
```

```
  // Runtime settings
```

```
runtime
```

```
  (Mode = 'Simulation', IsAccelerated = 'Yes', IsHybridExecution = 'No',  
   AuditOnSimulate = 'Yes', LogOnSimulate = 'Yes', EnableSimulationEmailNotification = 'No',  
   HibernateOnSimulate = 'No', EnableFixedStartTime = 'Yes',  
   SimulationStartTime = '1/28/2014 12:00 AM', AllowNewIterationsOnDeviceError = 'No',  
   EnableCongestionDetection = 'Yes', CongestionClearQueueTimeThreshold = '00:02:00',  
   MaxQueueTimeThreshold = '00:05:00', EnableVerboseLogging = 'No',  
   EnableExperiments = 'Yes', EnableCampaigns = 'Yes', AutoOffline = 'No',  
   SimplifiedInterface = 'No', ContainerLoadPrompting = 'Yes',  
   ContainerUnloadPrompting = 'Yes', AutoUnload = 'No') ;
```

```
// Devices and settings
```

```
devices
```

```
{
```

```
  KeyenceSR600 BarCode
```

```
    (OperationDefaults = '', Active = 'Active') ;
```

```
  MultidropCombi Combi
```

```
    (InitFluid = 'Default Fluid', PrimeWhenIdle = 'No',  
     PrimeOnInitialization = 'Yes', PrimeVolumeWhenIdle = '10',  
     PrimeIntervalWhenIdle = '5', CassetteUI = '1', OperationDefaults = '',  
     Active = 'Active') ;
```

```
  ContainerDataDriver ContainerDataDriver
```

```
    (ExcelPass = 'Momentum1', SummaryFormat = 'CSV', SummaryFilename = '',  
     SummaryColumns = 'DateTime,Location', OperationDefaults = '',  
     Active = 'Active') ;
```

```
  CytomatHotel CytomatHotel
```

```
    (ShakeDuringIncubate = 'No', RPMT1 = '100', RPMT2 = '100',  
     ScanBarcodeOnGet = 'Yes', FAMModeEnabled = 'No', SearchMode = 'Entire Device',  
     HotelsOccupancyLabel = '<Click to Edit ...>', ContainersParticipationLabel = '<Click to Edit ...>',  
     OperationDefaults = '', Active = 'Active') ;
```

DataMiner DataMiner

```
(OperationDefaults = '', Active = 'Active') ;
```

FreedomEVO EVO

```
(ProtocolPath = '\\\\Tecan-hp8300\\C\\ProgramData\\Tecan\\EVOware\\database\\scripts',
```

```
ParkScript = 'ThermoTest_Safe_Left.esc', ExecuteParkMethod = 'Yes',
```

```
UserName = 'Admin', UserPass = 'admin1', OperationDefaults = 'Thermo.Automation.Devices.Instruments.Tecan.FreedomEVO.RunScriptParameter\\Result~\\'\\';
```

```
Active = 'Active') ;
```

FileManager FileManager

```
(OperationDefaults = '', Active = 'Active') ;
```

FreeNest FreeNest

```
(OperationDefaults = '', Active = 'Inactive') ;
```

HiG4Centrifuge HiG4Centrifuge

```
(CounterWeight = 'Bucket2', OperationDefaults = '',
```

```
Active = 'Active') ;
```

Hotel Hotel_1

```
(OperationDefaults = '', Active = 'Active') ;
```

Hotel Hotel_2

```
(OperationDefaults = '', Active = 'Active') ;
```

InfiniteReader Infinite

```
(ProtocolPathListUI = '<Click Button to Edit>', KeepDoorOpenBetweenRead = 'No',
```

```
OperationDefaults = '', Active = 'Active') ;
```

IncuShake Inheco

```
(XFreq = '142', YFreq = '142', XAmp = '20', YAmp = '20',
```

```
XYPhaseShift = '0', OnlineTemp = '-1', OfflineTemp = '-1',
```

```
ShakeDuringLoadIncubate = 'No', UseExactDuration = 'No',
```

```
OperationDefaults = '', Active = 'Active') ;
```

LC480 LC480

```
(OperationDefaults = '', Active = 'Active') ;
```

MomentumOperator MomentumOperator

```
(OperationDefaults = '', Active = 'Active') ;
```

GenericMover Orbitor

```
(ParkLocation = 'STDloc:safe', ParkMoverAtEndOfRun = 'Yes',
```

```
MotionSettings = 'Velocity: 100%, Acceleration: 100%, Jerk: 100%',
```

```
AllowLidding = 'Yes', OperationDefaults = '', Active = 'Active') ;
```

PlateLoc PlateLoc

```
(PlateLocProfileName = 'Seal', PreHeatingTemperature = '160',
```

```
PostCoolOffTemperature = '40', CloseStageOnLoad = 'Yes',
```

```

    OperationDefaults = '', Active = 'Active') ;
Regrip Regrip
    (OperationDefaults = '', Active = 'Active') ;
AgilentMicroplateLabeler Vcode
    (ProfileName = 'LOR3413', OperationDefaults = '', Active = 'Active') ;
Waste Waste
    (OperationDefaults = '', Active = 'Active') ;
}

// Device Pools

pools
{
    StoragePool Hotels
        (Strategy = 'Priority First Available', ResetOnNewWorkUnit = 'Yes',
        SkipError = 'No', SkipOffline = 'No', OfflineSkipDuration = '00:00:00') Hotel_1, Hotel_2 ;
}

// Profile variables

variables
{
    Boolean Lock
        (DefaultValue = 'No', PromptForValue = 'No', Persist = 'No',
        Comments = '') ;
}

// ***** Version 1 *****

// User: Admin Date: Tuesday, June 03, 2014 Time: 10:31:20 PM

process [Read blank assay plates]
{
    // Containers

    containers
    {
        plate AssayPlate

```

```

    (WithLidOffset = '-5', MoverLiddingGripOffset = '3',
    WithLidHeight = '17', Thickness = '1', SealThickness = '0',
    Lid = '(None)', NumberOfWellRows = '8', NumberOfWellColumns = '12',
    WellNumberingMethod = 'Rows', ContainerTypeNameId = '4titude 4ti-0223',
    BarCodeRegularExpression = '', BarCodeFile = '',
    BarCodeAutoExpression = '"NC" + Format(Now, "yyMMddHHmmss") + "." + Format(WallClock, "fff")',
    GripOffset = 'Identity', GripForce = '0', Height = '15',
    StackHeight = '13.13', SetSize = '1', Attributes = '' ) ;
}

// Process steps

Infinite [Run Protocol]
    (ProtocolName = 'FluorescenceSingleCorning3561.mdfx',
    ProtocolPath = '\\\\Tecan-hp8300\\c\\Momentum Protocols\\Infinite',
    OutputNameFormat = '<PROTOCOL>_<DATETIME>', OutputPath = '\\\\Tecan-hp8300\\c\\Momentum Protocols\\Infinite',
    MaximumOperationTime = '00:30:00', PreCondition = '',
    RunOnAbortedIteration = 'No', Duration = '00:01:00',
    Comments = '', MinDelay = '00:00:00', MaxDelaySpecified = 'No',
    RequestedMaxDelay = '00:00:00', SpoilIfMaxDelayExceeded = 'No',
    Enabled = 'Yes', Result = '')
    AssayPlate 'Unlidded' in 'Nest' ;
}
}

```